



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Type K Retractable Curly Thermocouple Lead (ANSI)



XE-7031-001

TYPE-K (ANSI)



XE-7031-001

**TYPE-K
(ANSI)**

XE-7031-001

TYPE-K (ANSI)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Disponible en type K

Couleur codée a l'ANSI

2 metres de longueur lorsqu'il est entierement étendu

Terminé avec une prise miniature ajustée

Polyuréthane externe

Spécifications

Specifications

Product Code	XE-7031-001
Cable Length	2M fully extended
Cable Type	Polyurethane outer
Connector Type	Thermocouple Plug
Connector Size	Miniature
Colour	Yellow
Standards Met	Colour code (cores & jacket) to ANSI MC96.1
Thermocouple Type	K
Termination Type	Miniature Thermocouple Plug