




**LABFACILITY**

TEMPERATURE & PROCESS TECHNOLOGY

## Type K Retractable Curly Thermocouple Lead (ANSI)



XE-7031-001

**TYPE-K (ANSI)**



XE-7031-001

**TYPE-K (ANSI)**

XE-7031-001

**TYPE-K (ANSI)**

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Disponible en type K

Couleur codée a l'ANSI

2 metres de longueur lorsqu'il est entierement étendu

Terminé avec une prise miniature ajustée

Polyuréthane externe

### Spécifications

## Specifications

<b>Product Code</b>	XE-7031-001
<b>Cable Length</b>	2M fully extended
<b>Cable Type</b>	Polyurethane outer
<b>Connector Type</b>	Thermocouple Plug
<b>Connector Size</b>	Miniature
<b>Colour</b>	Yellow
<b>Standards Met</b>	Colour code (cores & jacket) to ANSI MC96.1
<b>Thermocouple Type</b>	K
<b>Termination Type</b>	Miniature Thermocouple Plug