




Type T PVC Isolating Mylar Screened Thermocouple Cable / Fil (IEC)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Type T PVC Isolé Mylar Screened Thermocouple Câble / Fil (IEC)

Le câble de compensation ressemble a la caractéristique thermoélectrique du thermocouple pertinent sur une plage de température ambiante limitée, de 0 a 80 oC en général. Étant donné que ces câbles sont fabriqués a partir de matériaux a faible cout, des économies peuvent etre réalisées sur les installations de l'usine par rapport a l'exécution du véritable câble d'extension de thermocouple.

Spécifications
Specifications

General Description	Good general purpose insulation for 'light' environments. Waterproof and very flexible.
Cable Length	10m, 25m, 50m, 100m, 200m
Cable Type	Twisted cores & screened construction for reduced electrical noise
Core / Strands	7/0.2mm
Conductor Cross Section 7/0.2 (mm²)	0.219
Outer Colour	Brown
+ Colour	Brown
- Colour	White
Drain Wire	Yes
Insulation Material	PVC Insulated with a mylar screen
Loop Resistance 7/0.2 (Ohms per combined metre)	2.3
Max. Temperature	+105°C
Min. Temperature	-10°C



+ Conductor	Copper
- Conductor	Constantan
Screen	Yes
Screen Type	Mylar
Standards Met	Colour code (cores & jacket) to IEC-584-3
Thermocouple Type	T
Tolerance	Class 2