

Raccords compression en laiton - filetage parallèle (BSPP)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Pour adapter une large gamme de capteurs isolation minrale et dautres capteurs de type sonde proposs par Labfacility.

Les raccords compression en laiton avec filetage parallle BSPP (British Standard Pipe Parallel) sont largement utilis pour connecter en toute scurit des tuyaux et des tubes dans diverses applications.

Les filetages BSPP ont un diamtre constant, ce qui ncessite une rondelle dtanchit ou un joint torique pour crer un joint tanche.

Fabriqus en laiton, ces raccords offrent une bonne rsistance la corrosion, une bonne durabilit et une facilit dinstallation, ce qui les rend adapts la plomberie, au CVC et aux applications basse et moyenne pression dans les systmes deau, dair et de gaz.

Spifications

Specifications

General Description	Male in-line stainless steel compression type with nut & olive
Material	Brass
Thread	1/8"BSP or 1/4"BSP available
Probe Diameter	Suitable for both metric and imperial probe diameters from either 1.5mm, 3.0mm or 3/16"