


## Raccords compression en laiton - Filetage conique (BSPT)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Pour adapter une large gamme de capteurs isolation minrale et dautres capteurs de type sonde proposs par Labfacility.

Les raccords compression en laiton avec filetage BSPT (British Standard Pipe Taper) sont largement utiliss pour connecter des tuyaux et des tubes dans diverses applications.

Les filetages BSPT sont coniques, ce qui signifie que le diamtre du filetage diminue vers la fin, ce qui amliore ltanchit lorsque le raccord est serr.

Fabriqus en laiton, ces raccords offrent une bonne rsistance la corrosion, une bonne durabilit et une facilit dinstallation, ce qui les rend adapts la plomberie, au CVC et aux systmes basse et moyenne pression impliquant de leau, de lair et du gaz.

### Spfications

## Specifications

<b>General Description</b>	Male in-line stainless steel compression type with nut & olive
<b>Material</b>	Brass
<b>Thread</b>	1/8"BSPT, 1/4"BSPT, 3/8"BSPT or 1/2"BSPT available
<b>Probe Diameter</b>	Suitable for both metric and imperial probe diameters from 1.5mm