

Casquettes de baÑfÂ-onnette - Acier inoxydable



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

General Description	Bayonet caps provide a method of quick fitting into suitable adapters located in the process; this technique is widely used in plastics machinery
Material	Stainless steel
Dimensions	11.5mm ID to 12.5mm ID
Sizes	to suit spring 3/16" or 6.0mm