



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Casquettes De BañfÂ-tonnette - BNP Brass



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

General Description	Bayonet caps provide a method of quick fitting into suitable adapters located in the process; this technique is widely used in plastics machinery
Material	BNP Brass
Dimensions	11.5mm ID to suit spring 3/16"