



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Adaptateurs de BaïfÂ-onnette - BNP Brass



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

General Description	Bayonet adaptors available in various lengths and thread sizes
Material	BNP Brass
Dimensions	1/8" BSPP and M10
Length	22, 35, 48, 64, 100 mm lengths