

## Adaptateurs de BañfÂ-onnette - BNP Brass



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Specifications

<b>General Description</b>	Bayonet adaptors available in various lengths and thread sizes
<b>Material</b>	BNP Brass
<b>Dimensions</b>	1/8" BSPP and M10
<b>Length</b>	22, 35, 48, 64, 100 mm lengths