


# Thermopockets haute pression percÃfÂ©s solides



- NOTES:  
 1. Solid Bar Stock Drilled Thermowell  
 2. Material 316/316L St.Stl  
 3. Pocket Rated 5000PSI  
 4. Material Certificate Required for each Batch

Dimensions in mm unless specified otherwise



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Specifications

<b>General Description</b>	Solid Drilled High Pressure Thermopockets
<b>Material</b>	316/316L Stainless steel
<b>Pocket Rated</b>	5000 PSI
<b>Dimensions</b>	316SS Tapered (12mm at tip)
<b>Length</b>	100mm, 150mm, 200mm stem length
<b>Type &amp; Stem Diameter</b>	316SS Tapered (12mm at tip)
<b>Stem Length</b>	100mm, 150mm, 250mm
<b>Standards Met</b>	Material Certificate supplied for each order