



## Thermopockets fabriquÃ©s



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Specifications

<b>General Description</b>	Fabricated Thermopockets in varying sizes for use with Industrial Style Temperature Probes.
<b>Material</b>	Stainless steel
<b>Length</b>	316SS with 7/16" or 8mm stem
<b>Type &amp; Stem Diameter</b>	316SS with 7/16" or 8mm stem
<b>Stem Length</b>	50mm, 100mm, 150mm, 350mm, 550mm
<b>Probe Diameter</b>	100mm, 150mm, 350mm, 550mm