




Thermocouple jonction expose fil nu de calibre fin (conducteurs de 0,025 mm) - Type K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Thermocouple ultra-fin pointe soude de type K pour une sensibilit optimale.

Spfications

Specifications

Product Code	XE-3517-001
General Description	Very fast response to a fraction of a second, Ideal for minimal invasive temperature sensing.
Sensor Type	Exposed Junction
Thermocouple Type	K
Cable Length	200mm
Cable Type	Fine gauge bare wire
Core / Strands	0.001" (0.025mm)
Max. Temperature	315°C