




# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## Tampons de fixation Thermo en tissu de verre PTFE



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Tampons auto-adhsifs en tissu de verre PTFE pour localiser les thermocouples pointe soude

Nominal -73C 260C

Taille du pad 13x18mm

La feuille contient 20 tampons

## Specifications

<b>Product Code</b>	XE-3599-001
<b>Series</b>	ATTACHMENT PADS