



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Plomb en caoutchouc de silicone Pt100 usage gnral - Type PRT - Diamtre 6,0 mm - Longueur 50 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-3630-001
General Description	Suitable for general purpose industrial & commercial use
Sensor Type	Pt100 (100 Ohms @ 0°C) to IEC 751, Class B, 4 wire
Probe Diameter	6.0mm
Probe Length	50mm
Sheath Material	316 stainless steel sheath
Cable Length	2M
Cable Type	Flexible silicone rubber insulated
Core / Strands	7/0.2mm
Cable Termination	Bare tails
Termination Type	4-wire tails allows connection to any PT100 instrument 2, 3 or 4 wire
Max. Temperature	200°C
Min. Temperature	-50°C
Accuracy	Class B
Number of Wires	4
Series	General Purpose Pt100