




# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## Thermomtre rsistance Pt100 4 fils de classe B, tte KNE avec longueur de sonde de 400 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Vous avez besoin de calibrer votre capteur de température ?**

Traable selon les normes dtalonnage ISO17025 (UKAS), Labfacility offre un choix dtalonnage traable en 3 ou 5 points laide de notre installation dtalonnage interne, vous pouvez slectionner nimporte quel point de température entre -10C et +200C.

**Veillez noter que si vous ajoutez une option de calibrage votre commande, vous serez contact pour les informations de livraison.**

**Spfications****Specifications**

<b>Product Code</b>	XE-3654-001
<b>Sensor Type</b>	Pt100 (100 Ohms @ 0°C) to IEC 751, Class B, 4 wire
<b>Probe Diameter</b>	6.0mm
<b>Probe Length</b>	400mm
<b>Head Type</b>	KNE
<b>Termination Type</b>	4 wire connection block, M20 x 1.5mm cable entry (gland included)
<b>Max. Temperature</b>	+450°C (connection head @ 170°C)
<b>Min. Temperature</b>	-50°C
<b>IP Rating</b>	IP67
<b>Accuracy</b>	Class B
<b>Number of Wires</b>	4