

## Capteur patch autocollant PT100 - Type PRT - Longueur 5M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Specifications****Specifications**

<b>Product Code</b>	XE-3661-001
<b>General Description</b>	Suitable for general purpose industrial & commercial use
<b>Sensor Type</b>	Silicone rubber patch with self adhesive foil backing
<b>Element</b>	Pt100
<b>Patch dimensions</b>	L40 x W15mm
<b>Cable Length</b>	5M
<b>Cable Type</b>	PTFE insulated twisted leads
<b>Core / Strands</b>	7/0.2mm
<b>Cable Termination</b>	Bare tails
<b>Termination Type</b>	4-wire tails allows connection to any PT100 instrument 2, 3 or 4 wire
<b>Max. Temperature</b>	150°C
<b>Min. Temperature</b>	-50°C
<b>Accuracy</b>	Class B
<b>Number of Wires</b>	4
<b>Series</b>	Patch Pt100