


## Type K IEC Angle droit - Pointe de disque en cuivre



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Gaine en acier inoxydable avec poigne moule et 2 metres de rallonge enroule termine par une fiche thermocouple miniature.

## Specifications

<b>Product Code</b>	XE-3817-001
<b>General Description</b>	Surface temperature spring loaded thermocouple with right angled probe and copper disc tip
<b>Thermocouple Type</b>	K
<b>Calibration</b>	IEC
<b>Cable Length</b>	2m
<b>Cable Termination</b>	Fitted Plug
<b>Max. Temperature</b>	+600°C
<b>Dimensions</b>	Probe O8mm(17mmL) Tip x 60mm Hot Leg and O6mm x 180mm Cold Leg with O4.5mm copper disc.