

## Sonde de temprature de surface - Type K IEC - 90 mm de long avec pointe de 10 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Gaine en acier inoxydable avec poigne moule et 2 metres de rallonge enroule termine par une fiche thermocouple miniature.

## Specifications

<b>Product Code</b>	XE-3818-001
<b>General Description</b>	Surface temperature fast response thermocouple with ceramic tip and coiled element
<b>Thermocouple Type</b>	K
<b>Calibration</b>	IEC
<b>Cable Length</b>	2m
<b>Cable Termination</b>	Fitted Plug
<b>Max. Temperature</b>	+900°C
<b>Dimensions</b>	Probe O6 x 90mm long stem with O12 x 10mm tip. Coiled element diameter 7mm
<b>Series</b>	Handheld Probe