



Thermocouple disque fil fin CEI, PTFE - Type K

XE-3515-001
**FINE WIRE DISC THERMOCOUPLE
TYPE K**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Thermocouple disque de cuivre avec bandes adhesives kapton haute temprature (13x18mm x 5) pour une fixation rapide et facile sur la surface/le processus.

Specifications

General Description	Fast response – ideal for surface measurement
Sensor Type	Copper disc, 6.35mm diameter x 0.3mm thick
Thermocouple Type	K
Cable Length	1m
Cable Type	PTFE insulated twisted wire
Core / Strands	1/0.2mm
Termination Type	Bare tails
Max. Temperature	+200°C
Min. Temperature	-0°C
Standards Met	IEC-584-3