


## Thermocouple jonction expose ANSI isol en fibre de verre avec prise miniature - Types K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Vous avez besoin de calibrer votre capteur de température ?**

Traable selon les normes dtalonnage ISO17025 (UKAS), Labfacility offre un choix dtalonnage traable en 3 ou 5 points laide de notre installation dtalonnage interne.

**Veillez noter que si vous ajoutez une option de calibrage votre commande, vous serez contact pour les informations de livraison.**

**Spifications**

**Specifications**

<b>General Description</b>	Fast response, welded exposed junction
<b>Sensor Type</b>	Good temperature resistance (not suitable in fluids)
<b>Thermocouple Type</b>	K
<b>Thermocouple Junction</b>	Exposed
<b>Cable Length</b>	300mm, 1, 2, 5 or 10 metre
<b>Cable Type</b>	glassfibre insulated flat pair cable
<b>Core / Strands</b>	1/0.3mm conductors
<b>Cable Termination</b>	Miniature Plug
<b>Max. Temperature</b>	+350°C (short periods up to 400°C)
<b>Min. Temperature</b>	-60°C
<b>Tolerance</b>	ANSI MC96.1 Class 2
<b>Standards Met</b>	ANSI MC96.1