

Capteurs de classe B Pt100 avec plomb isolé au téflon - Type PRT



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

Sensor Type	Pt100 element in a 3.0mm or 4.0mm diameter 316 stainless steel probe
Element	Pt100 element to IEC 751 Class B
Probe Diameter	3.0mm or 4.0mm
Probe Length	from 25mm to 125mm
Sheath Material	316 stainless steel probe
Cable Length	1m
Cable Type	Teflon® insulated twisted lead
Core / Strands	7/0.16mm
Cable Termination	Bare Tails
Max. Temperature	+200°C
Min. Temperature	-50°C
Standards Met	IEC 751
Accuracy	Class B