




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

100 DÃfÂ©tecteurs de blessures de fil ÃfÂ double ÃfÂ©lÃfÂ©ment Ohm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

General Description	Wirewound construction for optimum performance and stability
Sensor Type	Dual element Pt100
Max. Temperature	+650°C
Min. Temperature	-200°C
Tolerance	Class B or A to IEC 60751