




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

100 DÃfÂ©tecteurs de blessures de fil ÃfÂ double AfÂ©IAfÂ©ment Ohm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

| | |
|----------------------------|--|
| General Description | Wirewound construction for optimum performance and stability |
| Sensor Type | Dual element Pt100 |
| Max. Temperature | +650°C |
| Min. Temperature | -200°C |
| Tolerance | Class B or A to IEC 60751 |