

RÃfÂ©chauffeurs de buses de mica



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Ces produits sont disponibles jusqu'à épuisement des stocks avec des quantités limitées disponibles pour chaque version, vous serez alerté lors de l'ajout au panier si nous n'avons pas de stock suffisant -
[veuillez nous contacter](#) si vous avez des exigences spécifiques.

Specifications

General Description	Best quality Mica for high dielectric strength and maximum heat transfer
Element	High performance precision wound element for superior heat distribution
Cable Length	1 metre long stainless steel braided cable
Max. Temperature	Recommended for temperatures up to 300°C
Output	240V
Features #1	Full width low expansion straps
Features #2	Specially treated finish to resist oxidation and enhance thermal conductivity
Important Information	Wide range of applications including plastic process machinery, injection / blow moulding, rubber moulding, test equipment, food process industry etc