



Sonde thermocouple usage gnral CEI, cble isol en fibre de verre avec sur-tresse en acier inoxydable se terminant par des queues nues, fiche miniature ou standard - Type K, J



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

| | |
|------------------------------|--|
| General Description | Rugged thermocouple with stainless steel sheath |
| Thermocouple Type | K or J |
| Thermocouple Junction | Grounded hot junction for fast response |
| Probe Diameter | 4.5mm |
| Probe Length | 150mm |
| Cable Length | 1 or 2m |
| Cable Type | Glassfibre insulated lead with stainless steel overbraid |
| Core / Strands | 7/0.2mm |
| Termination Type | Bare tails, Miniature or Standard Plug (with cable clamps) |
| Max. Temperature | +350°C |
| Min. Temperature | -60°C |