

# Thermocouple isolante minérale de modèle industriel avec tête de terminal de KNE - Types K,J



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Besoin de votre capteur de température calibré?**

Traçable aux normes d'étalonnage ISO17025 (UKAS), Labfacility offre un choix d'un étalonnage traçable de 3 ou 5 points à l'aide de notre installation d'étalonnage interne, vous pouvez sélectionner tous les points de température entre -10 °C et 200 °C.

**Veillez noter que si vous ajoutez une option d'étalonnage à votre commande, vous serez contacté concernant les informations de livraison.**

## Specifications

<b>General Description</b>	Highly flexible, sheath can be bent/formed to suit many applications and processes
<b>Sensor Type</b>	Mineral insulated Thermocouple fitted with KNE terminal head
<b>Thermocouple Type</b>	Type K and J available
<b>Thermocouple Junction</b>	Insulated hot junction
<b>Probe Diameter</b>	6.0mm diameter
<b>Sheath Material</b>	310 stainless steel sheath
<b>Head Type</b>	KNE, Cast aluminium finish, screw type lid with chain
<b>Gland</b>	M16 x 1.5mm Plated brass gland fitted
<b>Max. Temperature</b>	+1100°C
<b>Min. Temperature</b>	-40°C
<b>IP Rating</b>	IP67