





# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## Queues nues isoles PFA ANSI - Type K,T



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Vous avez besoin de calibrer votre capteur de température ?**

Traable selon les normes dtalonnage ISO17025 (UKAS), Labfacility offre un choix dtalonnage traable en 3 ou 5 points laide de notre installation dtalonnage interne.

**Veillez noter que si vous ajoutez une option de calibrage votre commande, vous serez contact pour les informations de livraison.**

**Spifications****Specifications**

<b>General Description</b>	Fast response, welded exposed junction
<b>Thermocouple Type</b>	K or T
<b>Thermocouple Junction</b>	Exposed
<b>Cable Length</b>	300mm, 1, 2, 5 or 10 metre
<b>Cable Type</b>	Teflon® PFA insulated twin-twisted lead
<b>Core / Strands</b>	1/0.2, 1/0.3mm solid or 7/0.2mm stranded
<b>Cable Termination</b>	Bare wire
<b>Max. Temperature</b>	+250°C
<b>Min. Temperature</b>	-75°C
<b>Standards Met</b>	ANSI MC96.1