



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Certificat de conformitã©

ISO9001

**CERTIFICATE OF
CONFORMITY**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Un document certifié par une autorité compétente que le bien ou le service fourni répond aux spécifications requises.

Aussi appelé certificat de conformité ou certificat de conformité.

Specifications

Product Code	PO-COC
--------------	--------