


# Thermomtre rsistance hyginique Pt100, raccord de style RJT 1,5" avec metteur intgr - 3 fils, 6,0 mm x 75 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Thermomtre rsistance hyginique Pt100, raccord de style RJT 1,5"

## Specifications

<b>Product Code</b>	XE-5525-001
<b>General Description</b>	Ideal for hygienic processes including dairy, brewery and pharmaceutical applications.
<b>Sensor Type</b>	Pt100, 3 wire single element class A
<b>Head Type</b>	KNEF
<b>Connection Type</b>	Pt100, 3 wire single element class A
<b>Transmitter</b>	SEM-203p transmitter, ranged -50°C to +150°C
<b>Support Dia (D2)</b>	8.0mm
<b>Support Length (L2)</b>	50mm