



# Thermomtre rsistance hyginique Pt100, raccord de style RJT 1,5" sans metteur - 4 fils, 6,0 mm x 125 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Les raccords RJT (Ring Joint Type) sont souvent appels raccords standard. Le compromis dun joint torique qui laisse une petite crevasse o la doublure et la partie mle de lunion se chevauchent.

Ideal pour les processus hyginiques, y compris les applications laitires, brassicoles et pharmaceutiques.

### Spcifications

## Specifications

<b>Product Code</b>	XE-5524-001
<b>General Description</b>	Ideal for hygienic processes including dairy, brewery and pharmaceutical applications.
<b>Sensor Type</b>	Pt100, 4 wire single element class A
<b>Head Type</b>	KNEF
<b>Connection Type</b>	Pt100, 4 wire single element class A
<b>Transmitter</b>	None
<b>Support Dia (D2)</b>	8.0mm
<b>Support Length (L2)</b>	50mm