


Thermomtre rsistance hyginique Pt100, raccord de style RJT 1,5 " avec metteur intgr - 3 fils, 6,0 mm x 125 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Les raccords RJT (Ring Joint Type) sont souvent appels raccords standard. Le compromis dun joint torique qui laisse une petite crevasse o la doublure et la partie mle de lunion se chevauchent.

Idal pour les processus hyginiques, y compris les applications laitires, brassicoles et pharmaceutiques.

Spifications

Specifications

Product Code	XE-5526-001
General Description	Ideal for hygienic processes including dairy, brewery and pharmaceutical applications.
Sensor Type	Pt100, 3 wire single element class A
Head Type	KNEF
Connection Type	Pt100, 3 wire single element class A
Transmitter	SEM-203p transmitter, ranged -50°C to +150°C
Support Dia (D2)	8.0mm
Support Length (L2)	50mm