

## Thermomètre de résistance hygiénique Pt100, 1,5" Tri-Clamp



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Idéal pour les processus d'hygiène, y compris les applications laitières, brasseries et pharmaceutiques.

## Specifications

<b>Sensor Type</b>	Pt100, 3 or 4 wire single element class A
<b>Probe Diameter</b>	6mm
<b>Probe Length</b>	75mm OR 125mm
<b>Head Type</b>	KNEF
<b>Transmitter</b>	Optional SEM-203p transmitter, ranged -50°C to +150°C
<b>Number of Wires</b>	3