

Thermomtre rsistance hyginique Pt100, raccord 1,5" Tri-Clamp sans metteur - 4 fils, 6,0 mm x 75 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Idéal pour les processus hyginiques, y compris les applications laitieres, brassicoles et pharmaceutiques.

Spfications

Specifications

Product Code	XE-5527-001
Sensor Type	Pt100, 3 or 4 wire single element class A
Probe Diameter	6.0mm
Probe Length	75mm
Head Type	KNEF
Transmitter	None
Number of Wires	3