


**Thermomtre rsistance hyginique Pt100, raccord  
Tri-Clamp 1,5" avec metteur intgr - 3 fils, 6,0 mm x  
125 mm**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Idéal pour les processus hyginiques, y compris les applications laitieres, brassicoles et pharmaceutiques.

### Spfications

## Specifications

<b>Product Code</b>	XE-5530-001
<b>Sensor Type</b>	Pt100, 3 wire single element class A
<b>Probe Diameter</b>	6.0mm
<b>Probe Length</b>	125mm
<b>Head Type</b>	KNEF
<b>Transmitter</b>	SEM-203p transmitter, ranged -50°C to +150°C
<b>Number of Wires</b>	3