


## Thermocouple jonction expose CEI monoshot en PTFE avec fiche miniature - Type K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Un thermocouple isol en PTFE rponse trs rapide, fabriqu partir dune construction de cble unique injection unique qui fournit une isolation tanche au gaz, la vapeur et leau le long du cble.

Ce capteur est parfaitement adapt la surveillance de la temprature des pices dessai dans les applications dautoclave et aérospatiales.

| Besoin de calibrer votre capteur de temprature ?   |
|--|
| <p>Traceable aux talons dtalonnage ISO17025 (UKAS), Labfacility offre un choix dtalonnage traable en 3 ou 5 points laide de notre installation dtalonnage interne, vous pouvez slectionner nimporte quel point de temprature entre -10C et ; +200C.</p> <p><b>Veuillez noter que si vous ajoutez une option de calibrage votre commande, vous serez contact pour les informations de livraison</b></p> |

## Specifications

|                            |  |
|----------------------------|--|
| <b>General Description</b> | Unique single shot cable construction Resistant to chemicals, oils and fluids. Suitable for aerospace applications and is widely used in autoclaves. |
| <b>Sensor Type</b>         | Exposed Junction   |
| <b>Thermocouple Type</b>   | K & T  |
| <b>Cable Length</b>        | 300mm, 1, 2, 5 & 10m   |
| <b>Cable Type</b>          | PTFE   |
| <b>Core / Strands</b>      | 1/0.376mm  |
| <b>Termination Type</b>    | Terminated with a miniature flat pin plug (Plug max temp 220°C)  |
| <b>Max. Temperature</b>    | +250°C   |
| <b>Standards Met</b>       | IEC-584-3  |
| <b>Connector Type</b>      | Plug   |
| <b>Connector Size</b>      | Miniature  |