


PFA Cable Tidy IEC Exposed Junction Thermocouple avec mini-fiche intgre - Types K, T



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Thermocouple de rangement par cble avec prise IEC miniature installe, PFA double torsade ou paire plate.

Un thermocouple jonction expose fil fin de precision complet, haute sensibilit, avec un connecteur miniature intgr.

Vous avez besoin dtalonner votre capteur de temprature ?

Traable aux talons dtalonnage ISO17025 (UKAS), Labfacility offre le choix dun talonnage traable en 3 ou 5 points laide de notre installation dtalonnage interne, vous pouvez slectionner nimporte quel point de temprature entre -10 C et + 200 C.

Veillez noter que si vous ajoutez une option dtalonnage votre commande, vous serez contact concernant les informations de livraison.

Specifications

General Description	A complete, precision fine wire, high sensitivity thermocouple with a fitted miniature connector and cable tidy.
Sensor Type	Exposed Junction
Cable Length	1, 2m
Cable Type	PFA insulated wire, flat pair or twin twisted available
Termination Type	Miniature flat pin plug, colour coded in accordance IEC 584
Max. Temperature	Rated 220°C continuous
Tolerance	Class 1
Standards Met	IEC 584
Connector Type	Plug
Connector Size	Miniature