



Thermocouple jonction expose ANSI avec mini-prise intgre - Types K,T



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Vous avez besoin de calibrer votre capteur de température ?

Traable selon les normes dtalonnage ISO17025 (UKAS), Labfacility offre un choix dtalonnage traable en 3 ou 5 points laide de notre installation dtalonnage interne.

Veillez noter que si vous ajoutez une option de calibrage votre commande, vous serez contact pour les informations de livraison.

Spifications**Specifications**

General Description	A complete, precision fine wire, high sensitivity thermocouple with a fitted miniature connector
Sensor Type	Cable Tidy Thermocouple with Fitted Miniature ANSI Plug
Thermocouple Type	K or T
Thermocouple Junction	Exposed
Cable Length	1 or 2 metre
Cable Type	PFA twin twisted
Cable Termination	Miniature Plug
Max. Temperature	220°C
Tolerance	Class 1 wire for high accuracy
Standards Met	ANSI colour code