


## Thermocouple jonction expose ANSI avec mini-prise intgre - Types K,T



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Vous avez besoin de calibrer votre capteur de température ?**

Traable selon les normes dtalonnage ISO17025 (UKAS), Labfacility offre un choix dtalonnage traable en 3 ou 5 points laide de notre installation dtalonnage interne.

**Veillez noter que si vous ajoutez une option de calibrage votre commande, vous serez contact pour les informations de livraison.**

**Spifications****Specifications**

<b>General Description</b>	A complete, precision fine wire, high sensitivity thermocouple with a fitted miniature connector
<b>Sensor Type</b>	Cable Tidy Thermocouple with Fitted Miniature ANSI Plug
<b>Thermocouple Type</b>	K or T
<b>Thermocouple Junction</b>	Exposed
<b>Cable Length</b>	1 or 2 metre
<b>Cable Type</b>	PFA twin twisted
<b>Cable Termination</b>	Miniature Plug
<b>Max. Temperature</b>	220°C
<b>Tolerance</b>	Class 1 wire for high accuracy
<b>Standards Met</b>	ANSI colour code