

Raccords compression en acier inoxydable - filetage NPT (NPT)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Pour adapter une large gamme de capteurs isolation minrale et dautres capteurs de type sonde proposs par Labfacility.

Les raccords compression en acier inoxydable avec filetage NPT (National Pipe Thread) sont largement utilis pour connecter en toute scurit des tuyaux et des tubes dans diverses applications industrielles.

Les filetages NPT sont coniques, ce qui signifie que le diamtre du filetage diminue vers lextrmit, crant ainsi une tanchit plus tanche au fur et mesure que le raccord est serr.

Le matriau en acier inoxydable offre une excellente rsistance la corrosion, une durabilit et une rsistance, ce qui rend ces raccords adapts aux environnements difficiles, notamment le traitement chimique, le ptrole et le gaz, et les applications marines

Spfications

Specifications

General Description	Male in-line stainless steel compression type with nut & olive
Material	NPT - Stainless Steel Compression Fitting
Thread	1/8" NPT, 1/4" NPT, 1/2" NPT, 3/4"NPT or 3/8"NPT available
Probe Diameter	Suitable for both metric & imperial probe diameters from 1.0mm to 12.0mm