

Raccords compression en acier inoxydable - filetage NPT (NPT)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Pour adapter une large gamme de capteurs isolation minérale et d'autres capteurs de type sonde proposés par Labfacility.

Les raccords compression en acier inoxydable avec filetage NPT (National Pipe Thread) sont largement utilisés pour connecter en toute sécurité des tuyaux et des tubes dans diverses applications industrielles.

Les filetages NPT sont coniques, ce qui signifie que le diamètre du filetage diminue vers l'extrémité, créant ainsi une jointure plus étroite au fur et à mesure que le raccord est serré.

Le matériau en acier inoxydable offre une excellente résistance à la corrosion, une durabilité et une résistance, ce qui rend ces raccords adaptés aux environnements difficiles, notamment le traitement chimique, le pétrole et le gaz, et les applications marines.

Specifications

Specifications

General Description	Male in-line stainless steel compression type with nut & olive
Material	NPT - Stainless Steel Compression Fitting
Thread	1/8" NPT, 1/4" NPT, 1/2" NPT, 3/4"NPT or 3/8"NPT available
Probe Diameter	Suitable for both metric & imperial probe diameters from 1.0mm to 12.0mm