

Raccords compression en laiton - filetage NPT (NPT)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Pour adapter une large gamme de capteurs isolation minrale et dautres capteurs de type sonde proposs par Labfacility.

Les raccords compression en laiton avec filetage NPT (National Pipe Thread) sont couramment utiliss pour connecter des tuyaux et des tubes dans diverses applications.

Les filetages NPT sont coniques, ce qui cre une tanchit plus tanche au fur et mesure que le raccord est serr.

Le laiton, connu pour sa rsistance la corrosion, sa durabilit et son usinabilit, rend ces raccords adapts la plomberie, au CVC et aux applications basse et moyenne pression impliquant de leau, de lair et du gaz.

Spcifications

Specifications

General Description	Male in-line stainless steel compression type with nut & olive
Material	Brass
Thread	1/8" or 1/4" NPT available
Probe Diameter	To suite 3.0mm or 6.0mm probe diameters