





LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

4 Core Silicone Rubber Insulated, Section arrondie, câble de capteur PRT / Fil



25M



4 CORE
XE-2798-025

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

4 Core Silicone Rubber Insulated, Section arrondie, câble de capteur PRT / Fil

Peut être utilisé pour étendre jusqu'à 4 configurations RTD fil. Couramment utilisé avec les thermomètres Pt100, mais adapté avec d'autres types Pt130/500/1000 etc.

Spécifications**Specifications**

General Description	4 core silicone rubber insulated 7/0.2mm tin plated copper conductors, twisted together with silicone rubber outer jacket
Cable Length	10, 25m
Cable Type	Silicone rubber insulated
Core / Strands	7/0.2mm
Insulation Material	Silicone rubber insulated, tin plated copper conductors, twisted together with silicone rubber outer jacket
Max. Temperature	+200°C
Min. Temperature	-40°C
Number of Wires	4
Standards Met	Cores colour coded 2 x Red / 2 x White to IEC-751