

GM1311 Thermomtre alimentaire numrique professionnel haute performance



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

VEUILLEZ NOTER : Le dlai de livraison sera de 7 10 jours ouvrables lors du choix de l'etalonnage

GM1311 Thermomtre numrique professionnel haute performance pour aliments/barbecue/viande

Temperature ambiante : 0 40 C (32 104 F)

Humidit de lenvironnement : 10-90% HR

Tension nominale : 3 V

1 x thermomtre alimentaire

1 x manuel d'utilisation en anglais

1 x CR2032

Spcifications

Specifications

Product Code	XE-4002-001
Accuracy	-10 - 100°C: +/-1°C / 100 - 200°C: +/-1.5°C
Resolution	0.1°C / °F
Colour	White
Dimensions	9.29 in x 1.34 in x 0.71 in (23.6 cm x 3.4 cm x 1.8 cm). Display dimension:44x18mm
Function #1	LCD display
Function #2	Temperature lock
Function #3	Max/Min measuring
Function #4	Conversion between Celsius degree and Fahrenheit degree
Function #5	Auto power off
Max. Temperature	300°C (572°F)
Min. Temperature	-50°C (-58°F)
Probe Length	149.5 x 6mm



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Weight	Weight: 0.95 oz (27 g)
Power Supply	CR2032 button battery x 1 (included)