




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Thermocouple patch en silicone CEI, PFA avec prise miniature - Type K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

General Description	Resistant to many oils & chemicals
Sensor Type	Moulded silicone patch thermocouple, 30 x 15mm patch, self adhesive with aluminium foil backing.
Thermocouple Type	K
Cable Length	1 or 2 metre
Cable Type	PFA insulated twin-twist wire
Termination Type	Miniature flat pin plug, colour coded in accordance IEC 584
Max. Temperature	+150°C
Min. Temperature	-50°C