



## Thermocouple disque fil fin CEI, PTFE avec fiche miniature - Type K,J



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Bandes adhesives kapton haute temperature (13x18mm x 5) pour une fixation rapide et facile sur la surface/le processus. Disque en cuivre, 6,35 mm de diamtre x 0,3 mm dpaisseur.

## Specifications

<b>General Description</b>	Exposed junction thermocouple with Copper Disc brazed at tip
<b>Sensor Type</b>	Copper disc
<b>Thermocouple Type</b>	Type 'K' or 'J' IEC, Type 'K' BS
<b>Thermocouple Junction</b>	Grounded measuring junction
<b>Cable Length</b>	1 or 2 metre
<b>Cable Type</b>	1/0.2mm PTFE insulated twisted wire
<b>Termination Type</b>	Miniature flat pin plug, colour coded in accordance IEC 584