




## Thermocouple baonnette rglable CEI, cble surtress en acier inoxydable en fibre de verre - Type K,J



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Thermocouple baonnette rglable Type K ou J IEC & Type J BS**

Thermocouple baonnette rglable avec ressort de 200 mm avec capuchon de fixation rglable (11,5 mm de diamtre intrieur, fente unique).

**Spcifications**
**Specifications**

<b>General Description</b>	Thermocouple with adjustable insertion depth
<b>Sensor Type</b>	6.0mm diameter tip
<b>Thermocouple Type</b>	K, J IEC J BS
<b>Thermocouple Junction</b>	Grounded tip for fast response
<b>Cable Length</b>	2, 3m
<b>Cable Type</b>	Glassfibre insulated lead with stainless steel overbraid
<b>Core / Strands</b>	7/0.2mm
<b>Cable Termination</b>	Bare wire or standard plug
<b>Max. Temperature</b>	+350°C
<b>Min. Temperature</b>	-60°C
<b>Standards Met</b>	IEC-584-3 / BS
<b>Bayonet cap</b>	11.5mm I.D. single slot
<b>Connector Type</b>	Plug
<b>Connector Size</b>	Standard
<b>Dimensions</b>	25x13mm shim (0.2mm thick)



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

<b>Thread</b>	M6 x 1mm thread
<b>Spring</b>	200mm