

## Master Class Thermomtre de cuisson pliant



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**VEUILLEZ NOTER: Le dlai sera de 7 10 jours ouvrables lors du choix de l'etalonnage**

Parfait pour la cuisson prcise du bnf, du veau, de l'agneau, du porc, du poulet et de la dinde, assurez-vous d'avoir un dlicieux rti de viande pour votre djeuner du dimanche, impressionnant les amis et la famille.

L'cran LCD facile lire montre clairement la temperature interne mesuree par la sonde en acier inoxydable, avec une plage de temperature de -45 C 200 C et de -49 F 392 F.

Ce thermomtre pliant est galement idal pour ranger facilement dans les tiroirs de cuisine, le mecanisme de pliage rduisant la taille globale.

### Spfications

## Specifications

<b>Product Code</b>	XE-3966-001
<b>General Description</b>	The easy to read LCD display clearly shows the internal temperature measured by the stainless steel probe
<b>Sensor Type</b>	Folding Stainless Steel Probe
<b>Max. Temperature</b>	200°C (392°F)
<b>Min. Temperature</b>	-45°C (-49°F)
<b>Dimensions</b>	25.8 x 8.8 x 2.2 cm
<b>Power Supply</b>	1 x LR44 battery (included)
<b>Weight</b>	100g