

Master Class Thermomtre de cuisson pliant



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

VEUILLEZ NOTER: Le dlai sera de 7 10 jours ouvrables lors du choix de l'etalonnage

Parfait pour la cuisson precise du bœuf, du veau, de l'agneau, du porc, du poulet et de la dinde, assurez-vous d'avoir un délicieux rti de viande pour votre djeuner du dimanche, impressionnant les amis et la famille.

L'cran LCD facile lire montre clairement la temperature interne mesuree par la sonde en acier inoxydable, avec une plage de temperature de -45 C 200 C et de -49 F 392 F.

Ce thermomtre pliant est galement idal pour ranger facilement dans les tiroirs de cuisine, le mecanisme de pliage rduisant la taille globale.

Spifications

Specifications

Product Code	XE-3966-001
General Description	The easy to read LCD display clearly shows the internal temperature measured by the stainless steel probe
Sensor Type	Folding Stainless Steel Probe
Max. Temperature	200°C (392°F)
Min. Temperature	-45°C (-49°F)
Dimensions	25.8 x 8.8 x 2.2 cm
Power Supply	1 x LR44 battery (included)
Weight	100g