

GM816 Anemomtre Anmomtre Anmomtre / Anmomtre numrique



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

GM816 Anmomtre Anmomtre Anmomtre / Anmomtre numrique

Ideal pour les utilisations industrielles et domestiques telles que la mesure de la vitesse du vent et de la temprature du processeur, des ventilateurs dordinateur, des climatiseurs, des souffleurs dair, etc.

Specifications

Product Code	XE-4003-001
Model	GM816 Anemometer Wind Meter Anemometer / Digital Anemometer
Accuracy	+/- 5%
Air Temp: Accuracy	+/- 2°C, +/- 3.6F
Air Temp: Range	-10~45°C, 14~113°F
Air Temp: Resolution	Resolution 0.2°C, 0.36°F
Air Velocity: Range	0.3 - 30m/s, 0 - 90km/h, 0 - 5860ft/min, 0 - 65mph, 0 - 55Knots
Air Velocity: Resolution	0.1m/s, 0.3km/h, 19ft/min, 0.2mph, 0.2Knots
Air Velocity: Threshold	0.1m/s, 0.3km/h, 39ft/min, 0.2mph, 0.1Knots
Battery	1 x 3V CR2032 battery (Not Included)
Colour	Yellow / Black
Current Consumption	3mA Approx
Dimensions	10.5 x 4.5 x 2 cm / (4.13 x 1.77 x 0.79)
Operating Humidity	less than or equal to 90%RH
Operating Temperature	-10~45°C (14~113°F)
Weight	105g / 3.7oz
Storage Temperature	-40~60°C, -40~140°F
Thermometer	NTC Thermometer
Series	Anemometer