


# Prise de connecteur de thermocouple de montage sur circuit imprimé IM-K-PCB Type K IEC avec broches plaques tain



**Miniature PCB Socket**  
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Prise de connecteur de thermocouple de montage sur circuit imprim IM-K-PCB Type K IEC avec broches plaques tain**

Ces articles ne peuvent pas tre commands en ligne, veuillez donc appeler ou envoyer un e-mail notre quipe de vente avec la quantite requise pour obtenir des informations sur le prix et la livraison.

Spicalement conu pour les circuits imprimis, broches solides plaques tain pour le montage direct sur la carte et mise disposition pour capteur CJC

**Spcifications****Specifications**

<b>Product Code</b>	XE-1109-TIN
<b>General Description</b>	Versatile, high quality circuit board mounting miniature socket
<b>Thermocouple Type</b>	K IEC
<b>Connector Type</b>	PCB Socket
<b>Connector Size</b>	Miniature
<b>Colour</b>	Green
<b>Max. Temperature</b>	220°C
<b>Pins</b>	Tin Plated
<b>+Positive Contact</b>	Nickel Chromium
<b>-Negative Contact</b>	Nickel Alloy
<b>Standards Met</b>	BS EN 50212:1997. RoHS